



Invacare® | **Propad™ Original Cushion**

Cut foam with castellations for comfort and support



The *Invacare Propad Original Cushion* has been designed to facilitate good pressure redistribution when used with both wheelchairs and standard armchairs. This cushion is constructed to ensure improved support and prolonged longevity.



Designed to deliver enhanced pressure care

The *Invacare Propad Original Cushion* is made from a layer of castellated foam, each segment delivers enhanced support and pressure area care to the individual. The “keyhole” cut at the base of each castellation ensures good air circulation, reducing potential for skin maceration.

Seating surface

Offers improved stability and postural support; castellations reduce shear and friction forces on tissue.

Features and Options



Cover helps prevent a build up of shear and friction forces

Cover is two-way stretch, water-resistant and vapour-permeable. Cover features an anti-slip base.



Available in 116 sizes

The *Invacare Propad Original Cushion* is available in 116 sizes ranging from 10-24” width and 10-22” depths (See prescription chart for details).



Full length flap conceals zip

The full length flap concealing the zip reduces the potential for fluid ingress into the cushion and protects the foam core. The cover can be removed and laundered.

Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please visit www.invacare.co.uk

	Size range	Height	Total product weight	Max. user weight	Max. cleaning temperature	Warranty	
Propad Original Cushion	A range of 116 sizes	102 mm	1 kg	108 kg (17 stone)	80°	1 year	Manufactured to European Standards: Interior Foam: BS EN 597-1 Cig, BS EN597-2 Match, BS 5852 (Crib 5) Cover: BS 7175 (Crib 5)

Colours

Fuscia

Black Smooth

