



## THE KÜSCHALL CHAMPION: AN EXAMPLE OF THE PIONEERING DEVELOPMENT OF KÜSCHALL AND ITS RELATION TO THE DAILY CHALLENGES RAINER KÜSCHALL IS FACING IN HIS DAY TO DAY LIFE.

Behind every küschall wheelchair development there is a story! After his swimming accident at the age of 16, Rainer Küschall lost all mobility and freedom. Despite the fear, he was always determined to overcome the daily challenges he faced. One of these was to transfer from his wheelchair to his car. At that time, he always needed another person to put the wheelchair in and out of the car but at least he had the freedom to drive... until a warm summer day when he started to feel unwell. Due to his injury, Rainer cannot sweat, which creates a huge risk to his health. Realising the immediate danger, he pulled over at a petrol station, opened the car door, and waved at two people standing outside in a desperate attempt to get some water.

When the people saw him, they were so afraid of his disability, that they immediately closed the car door again. Rainer was trapped in his own car! That was the moment Rainer became determined to take the next step to independence; making a wheelchair that he would be able to fold, move over his lap and put on the passenger seat. Thus, the idea of the küschall Champion was born.

Twenty-eight years on, the küschall Champion is the only wheelchair model of its type, with its tried and trusted design that has stood the test of time making it as popular today as when it was first launched.






## KÜSCHALL CHAMPION, THIS UNIQUE FOLDABLE WHEELCHAIR, IS NOW EVEN BETTER!

New features are the high stability and stiffness. The design and driving performance of this wheelchair is without doubt comparable to that of a rigid one. With stepless adjustments at the front and the rear, and easy positioning of the centre of gravity, the küschall Champion suits all types of individuals and allows them to configure the wheelchair to their needs.

The new design and the smart construction and assembly solutions make this wheelchair feel like a rigid one.

Its unique horizontal folding mechanism turns this wheelchair into one of the smallest possible active chairs, unrivalled in the market since 1986. Only one movement is needed to fold the wheelchair and its backrest. Compact and lightweight, its extremely easy to fold and handle; ideal for active individuals who are constantly on the go, and it fits neatly into a small car.

Küschall Champion is available in Aluminium, Titanium and Carbon.





THE KÜSCHALL CHAMPION OFFERS A WIDE CHOICE OF COLOURS, CONFIGURATIONS, AND OPTIONS ALLOWING YOU TO PERSONALISE THIS WHEELCHAIR AS YOU WISH.

EXPRESS YOURSELF WITH A WHEELCHAIR THAT PERFECTLY FITS



YOUR PERSONALITY.









ARE YOU LOOKING FOR A WHEELCHAIR WITH STYLE AND DURABILITY? THE KÜSCHALL CHAMPION IS AVAILABLE WITH A TITANIUM OR CARBON FRAME. STRONG AND SCRATCH-RESISTANT, TITANIUM AND CARBON ARE THE PERFECT CHOICE FOR ACTIVE LIFESTYLES AND INDIVIDUALS WITHOUT ANY COMPROMISE ON DESIGN – A REAL EYE-CATCHER!







Colour Power Sets (CPS) are available in 6 different tones: chrome, technoanthracite, comet red, champagne, b-blue and arctic blue.



Compact size (765x390x285 mm): fold the wheelchair and its backrest, remove the rear wheels, and it fits in the smallest of cars.



There is a choice of 9 different seat positions and 5 backrest angles. Side supporters also come in CPS when selected.



As a standard feature it includes a soft locking mechanism on the backrest so that it doesn't open when folded.



Frame angle 75° / 85° Al & Ti (Straight or tapered) 85° Ca



360 - 460 mm Al, Ti & Ca 400 - 480 mm Tapered in increments of 20 mm



340 - 480 mm Al 340 - 460 mm Ti 400 - 460 mm Ca in increments of 20 mm



Seat height front: 450 – 540 mm rear: 390 – 490 mm in increments of 10 mm



Knee-to-heel length 320 – 500 mm in increments of 10 mm



Backrest height 300 – 465 mm in increments of 15 mm



**Total weight, unloaded** 9,2 kg



Transport weight (w/out rear wheels)



Max user weight 120 kg **Küschall AG** | Benkenstrasse 260 | CH-4108 Witterswil | Switzerland Phone +41 (0) 61 487 70 70 | Fax +41 (0) 61 487 70 79 | www.kuschall.com

**UNITED KINGDOM: Invacare Ltd.** | Pencoed Technology Park, Pencoed, Bridgend CF35 5AQ Phone +44 (0) 1656 776222 | Fax +44 (0) 1656 776220 | www.invacare.co.uk

**IRELAND: Invacare Ltd.** | Unit 5, Seatown Business Campus | Seatown Road | Swords, Co. Dublin Phone +353 (0)1 810 7084 | Fax +353 (0)1 810 7085 | www.invacare.ie

**Invacare EASTERN EUROPE, MIDDLE EAST, RUSSIA & CIS** | Kleiststraße 49 | D-32457 Porta Westfalica Phone +49 (0)57 31 754 540 | Fax +49 (0)57 31 754 541 | edo@invacare.com

Küschall is a registered trademark. Every effort has been made to ensure that the contents of this publication are fully up-to-date at the time of going to print. As part of its ongoing improvement of products, Küschall AG reserves the right to modify existing models at any time. Any use of this publication or of parts thereof as well as any reproduction of images must have the written consent of Küschall AG.

2015 Copyright © 01/2015, Küschall AG, Switzerland – All rights reserved

