



Invacare Birdie/Birdie Compact

Product Data Sheet

Product Range Available

Code: 3000507E2A00-0125 - Birdie Mobile Hoist, electrical legs, 2-point spreader bar, detachable battery

Code: 3000504D2A00-0125 - Birdie Mobile Hoist, manual legs, 2-point spreader bar, fixed battery

Code: 3000505D2A00-0125 - Birdie Mobile Hoist, manual legs, 2-point spreader, detachable battery

Code: 3000510D2A00-0125 - Birdie Compact Mobile Hoist, manual legs, 2-point spreader bar, fixed battery

Code: 3000513D2A00-0125 - Birdie Compact Mobile Hoist, manual legs, 2-point spreader bar detachable battery

All products are delivered with battery, charger lead and 2 point spreader bar as standard.

Detachable battery control boxes have Linak Jumbo Care™ functionality.

Key Features

- Folding format with magnetic locking pin
- No tools for assembly
- SWL - Birdie 180kg, Birdie Compact 150kg
- 360 degree range even in highest position
- Steel construction
- Externally tested and verified to ISO 10535 by TÜV
- Manual and Electrical leg versions
- Warranty 5 years*

Technical Specification

- 150kg (23.5 stone) and 180 kg (28 stone) weight limit
- Overall product weight - Birdie 42kg, Birdie Compact 36kg
- Castors 100mm (UK) (2 braked castors at rear)
- Maximum leg opening – Birdie 1040mm, Birdie Compact 890mm (internal)
- Closed leg dimension – Birdie 580mm, Birdie Compact 460mm (internal)
- Closed leg dimension – Birdie 655mm, Birdie Compact 530mm (external)
- Lift height range – Birdie 450 – 1670mm, Birdie Compact 500 – 1570mm
- Base length – Birdie 1150mm, Birdie Compact 1000mm
- Turning radius – Birdie 1400mm, Birdie Compact 1070mm
- Number of cycles per full charge – 40
- Battery – 2.9ah
- Emergency stop



- Emergency lower – Birdie Electrical and Manual, Birdie Compact Electrical
- **Accessories:**
- 2 point and 4 point spreader bars in 450 and 550mm widths
- Handle for leg spread
- Homecare scale
- Jumbo Care service interface diagnostic (for control boxes with Jumbo Care functionality)
- Wide range of slings
- Compatible with other manufacturers slings - for sling compatibility see bulletin PU661

Supporting Documentation

- Birdie range brochure
- Birdie User manual
- Weigh scale log book
- Jumbo Care installation manual
- Jumbo Care service data tool
- Sling and Hoist Assessment Guidelines

All above documentation can be found on www.invacare.co.uk

Invacare's Business Development Managers can provide in field support and training on hoists and slings

*5 years on control box, 3 years on frame and 1 year on battery and handset

Manufacturers Contact Details

Invacare Ltd, Unit 4 Pencoed Technology Park, Bridgend, CF35 5AQ

Tel: 01656 776222 **E-mail:** UK@invacare.com

Web: www.invacare.co.uk / www.thinkindependence.co.uk



| safe patient handling