# Invacare® | **Pronto™ M61**





# A fresh approach to mobility

The sleek and stylish *Invacare* **Pronto M61** is an affordable centre wheel drive power chair with 6 wheels on the ground for optimum safety and stability. It features *Invacare's* **Sure Step**\* technology for smoother driving over lower thresholds. Designed around a compact base, the **Pronto M61** provides excellent manoeuvrability especially in confined spaces, and the electric seat provides 12 cm of elevation.





#### **Optimum Comfort**


The padded, shaped seat ensures you are comfortable for day long use, and the simple to operate seat will electrically raise you up to 12 cm. In the elevated position the *Invacare Pronto M61* provides increased accessibility, which is ideal in office environment or when you are out and about. For easier transfers into or out of the **Pronto M61**, the seat will rotate 90 degrees in either direction, and the arm rest can flip up or be removed altogether. The foot plate is centre mounted and height adjustable to provide maximum comfort.

### Upwardly mobile

With its centre wheel drive and its compact dimensions, the **Pronto M61** gives excellent manoeuvring in very small spaces and has a very tight turning circle. Invacare's Sure **Step**\* technology gives smooth ride over low thresholds of up to 5 cm while maintaining optimum stability at all times. The powerful motors lift the front castors to assist changes in height for a smoother, less jarring ride over rough surfaces ensuring safety at all times.

## Taking it with you

Intuitive driving can be obtained with the true centre wheel drive system combined with the programmable Shark controls. This compact power chair can be easily dismantled for transportation in the boot of a car so you can take it with you wherever you go.



#### **Technical data**



**Pronto M61** 455 - 500 mm



410 - 460 mm



470 -600 mm 500 - 630 mm<sup>2</sup>





630 mm (excluding joystick)



790 mm 1050 mm 4

Invacare Limited



1020 mm



1000 mm







2 x 32 Ah



2 x 180 W



Pronto M61

Pronto M61

up to 19 km



6 km/h





136 ka



#### Frame colours



Deep blue pearl



Candy red pearl

More detailed information, technical data, references for transport or description of certain options can be found in the wheelchair's user manual.

- 1/ Seat height is measured including cushion, with 4º tilt angle
- 2/ Seat height is measured including cushion. with 9º tilt angle
- 3/ With footboard up
- 4/ With footboard down, with side mounted legrests
- 5/ Compliant with the directive on the Restriction of the use of certain Hazardous Substances in electrical and electronic equipment - RoHS (2011/65/EC).
- 6/ Compliant with the regulation concerning the Registration, Evaluation, Authorization and Restriction of Chemicals.



Follow us on YouTube: Invacare Ltd - Think Mobility

Pencoed Technology Park - Pencoed -Bridgend CF35 5AQ - United Kingdom Tel: +44 1 656 776 222 Fax: +44 1656 776 220 E-mail: uk@invacare.com Sales Order E-mail: ordersuk@invacare.com www.invacare.co.uk Invacare Ireland Ltd Unit 5 - Seatown Business Campus Seatown Road - SWORDS - County Dublin - Ireland Tel. +353 1 8107084 Fax +353 1 8107085 Email: Ireland@invacare.com www.invacare.ie © 2016 Invacare International Sàrl All rights reserved.

All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation.

ProntoM61 - UK - 01/2016 - 1468372





